

Universal & high speed 4 x 100 ml Z 323 /Z 323 K

Features

- Brushless drive
- Microprocessor control
- Imbalance switch off
- Lid lock
- Automatic Lid lock
- Automatic rotor recognition

Digital Display for:

- Speed
- Force
- Running time
- Deceleration
- Temperature
- Program memory
- Key for short runs
- Acceleration/deceleration profile
- Continuous run function
- Standstill cooling
- Preselection by potentiometers
- Diagnosis programme

Specifications Z 323

Max. volume:	4 x 100ml
Max speed:	17000 rpm
Max. RZB:	27464 xg
Dimensions	
w x h x d (cm):	39 x 37 x 55
Weight:	44.0 kg
Product no	12370/50

Specifications Z 323 K

Max. volume:	4 x 100ml
Max speed:	17000 rpm
Max. RCF:	27464 xg
Dimensions	
w x h x d (cm):	69 x 37 x 55
Weight:	79.0 kg
Product no	12370/51



Accessories

Product no	Description			
12370/52	Swing out rotor	4 x 100ml	4500 rpm	3350xg
12370/53		8 x 15ml Falcon		
12370/54		36 x 15ml		
12370/55	Swing out rotor	2 x 3 Microtitre plates	4500 rpm	2762 xg
12370/56	Angle rotor	6 x 50ml	6000 rpm	3904 xg
12370/57	Angle rotor	12 x 15ml	6000 rpm	4185 xg
12370/58	Angle rotor	24 x 15ml	3500 rpm	2000 xg
12370/59	Angle rotor	44 x 1.5/2.0ml	13500 rpm	17320 xg
12370/60	Hematocrit rotor	24 x capill. 1.3 x 75mm		14811 xg
12370/61	Angle rotor	6 x 85ml	13500 rpm	20987 xg
12370/62	Angle rotor	24 x 1.5/2.0ml	17000rpm	27464 xg

LAB