

Centrifuges

Universal 4 x 500 I Z 400 / Z 400 K



Features

- Brushless drive
- Microprocessor control
- Imbalance switch off
- Lid lock
- Automatic rotor recognition

Digital Display for:

- Speed
- Force
- Deceleration
- Temperature
- Key for short runs
- Acceleration/deceleration profile
- Continuous run function
- Standstill cooling
- Preselection by potentiometers
- Diagnosis programme

Specifications Z 400

Max. volume:	4 x 500ml
Max speed:	135000 rpm
Max. RZB:	17320 xg
Dimensions	
w x h x d (cm):	43 x 37 x 54
Weight:	39.0 kg
Product no	12370/70

Specifications Z 400 K

Refrigerated

Max. volume:	4 x 500ml
Max speed:	135000 rpm
Max. RZB:	17320 xg
Dimensions	
w x h x d (cm):	74 x 34 x 58
Weight:	89.0 kg
Product no	12370/71

Accessories

Product no	Description			
12370/72	Swing out rotor	4 x 500ml 48 x 15ml Falcon 28 x 50ml Falcon	3500 rpm	2547 xg
12370/73	Swing out rotor	2 x 3 Microtitre plates	4500 rpm	2762 xg
12370/74	Angle rotor	6 x 50ml	6000 rpm	3904 xg
12370/75	Angle rotor	12 x 15ml	6000 rpm	4185 xg
12370/76	Angle rotor	24 x 15ml	3500 rpm	2000 xg
12370/77	Angle rotor	44 x 1.5/.2.0ml	13500 rpm	17320 xg
12370/78	Hematocrit rotor	24 x capill. 1.3 x 75mm		14811 xg