

Universal & high speed 4 x 500ml Z 383 / Z 383 K



Features

- Brushless drive
- Microprocessor control
- Imbalance switch off
- Lid lock
- Automatic Lid lock
- Automatic rotor recognition

Digital Display for:

- Speed
- Force
- Running time
- Deceleration
- Temperature
- Program memory
- Key for short runs
- Continuous run function
- Standstill cooling
- Preselection by potentiometers
- Diagnosis programme

Specifications Z 383

Max. volume:	4 x 500ml
Max speed:	17000 rpm
Max. RZB:	27464 xg
Dimensions	
w x h x d (cm):	50 x 38 x 59
Weight:	69.0 kg
Product no	12370/80

Specifications Z 383 K

Refrigerated

Max. volume:	4 x 500ml
Max speed:	17000 rpm
Max. RCF:	27464 xg
Dimensions	
w x h x d (cm):	76 x 38 x 60
Weight:	107.0 kg
Product no	12370/81

Accessories

Product no	Description			
12370/82	Swing out rotor	4 x 500ml 48 x 15ml Falcon 28 x 50 ml Falcon 4 x 3 Microtitre plates	4500 rpm	3373xg
12370/83	Swing out rotor	2 x 3 Microtitre plates	4500 rpm	2762 xg
12370/84	Angle rotor	6 x 50ml	6000 rpm	3904 xg
12370/85	Angle rotor	12 x 15ml	6000 rpm	4185 xg
12370/86	Angle rotor	24 x 15ml	3500 rpm	2000 xg
12370/87	Angle rotor	44 x 1.5/.2.0ml	13500 rpm	17320 xg
12370/88	Angle rotor	24 x 1.5/.2.0ml	17000rpm	27464 xg
12370/89	Hematocrit rotor	24 x capill. 1.3 x 75mm		14811 xg
12370/90	Angle rotor	6 x 85ml	15000 rpm	25910 xg