

Flow filtration Sieva - 2



Specifications Sieva - 2

Max. volume:	500ml
Max speed:	10000 rpm
Max. RCF:	7825 xg
Dimensions	
w x h x d (cm):	39 x 36 x 49
Weight:	25.0 kg

Product no **12371/10**

Features

- Brushless drive
 - Microprocessor control
 - Imbalance switch off
 - Lid lock
- Digital Display for:**
- Speed
 - Force
 - Running time
 - Key for short runs
 - Acceleration/deceleration profile
 - Continuous run function
 - Preselection by potentiometers
 - Diagnosis programme

Accessories

Product no	Description				
12371/11	Centrifugal basket	500ml	perforated	10000 rpm	7825 xg
12371/12	Centrifugal basket	500ml	non-perforated	10000 rpm	7825 xg

Possibility for external temperature equalisation