

Centrifuge Accessories

Centrifuge Bottles - cont

Polycarbonate

Transparent and autoclavable, these bottles have wide mouths, and polypropylene sealing caps (with silicone gaskets) for excellent chemical resistance. Supplied in packs of four.



Product No	Capacity-ml	RCF (xG)
10859/12	250	27500
10859/13	450	13700

Polypropylene Copolymer

As 10859/12 & 13 bottles/caps above, but bottle material is translucent polypropylene copolymer. Supplied in packs of four.



Product No	Capacity-ml	RCF (xG)
10859/14	250	27500
10859/15	450	13700

Sealing Caps For above Centrifuge Bottles


Autoclavable caps, recommended to give extra security when spinning hazardous materials. 11056/01 are polypropylene closures with Tefzel ETFE plug, Viton fluoro-elastomer O-ring and shoulder support ring to reinforce bottle neck. 11056/02 are polypropylene closures with silicone gasket. Supplied in packs of two.



Product No	To Suit 10859/
11056/01	06, 07, 11
11056/02	04, 05

Centrifuge Tubes, Glass

A range of glass, rimmed centrifuge tubes. Supplied in packs of 10.



Product No	Size-ml	Type	Glass	Graduated
10864/01	15	A	Soda	No
10864/02	15	B	Soda	No
10864/03	15	C	Soda	Yes
10864/04	15	D	Borosilicate	No
10864/05	110x117mm	E	Soda	No

High quality impact and heat resistant borosilicate glass round bottom centrifuge tubes in the capacities and dimensions indicated below. These tubes conform to DIN 58970 where appropriate, and are rated for up to 4000 rpm.



Product No	OD-mm	Height-mm	Cap-ml	Pack Qty
10865/01	12	100	6	50
10865/02	16	100	12	50
10865/03	24	100	25	10
10865/04	34	100	50	10
10865/05	40	115	80	10
10865/06	44	100	80	10
10865/07	56	147	250	10

