

Laboratory Cleaning

DECON 90

A high quality cleaning agent which will remove stubborn contaminants such as agar, grease, blood, protein, radioactive deposits and many more, from glass, ceramics, plastics (NOT polycarbonate), stainless steel & ferrous metals. Decon 90 is non-toxic, phosphate-free and biodegradable containing no hypochlorite, E.D.T.A. or N.T.A. Economical in use, typically being diluted to 2-5%. NOT suitable for use on non-ferrous metals (eg, aluminium, zinc) for which Neutracon should be used



Product No	Description
10012/01	1 litre bottle
10012/02	5 litre bottle
10012/03	20 litre drum

Neutracon

An alternative cleaning agent for use on pH-sensitive materials (aluminium, copper, zinc, etc). Does NOT contain phosphates, enzymes, E.D.T.A., N.T.A. or chlorine



Product No	Description
10014/02	5 litre bottle

Virkon Disinfectant

Non-chlorine based disinfectant which is environmentally friendly and harmless in normal use. Effective against Aids, Hepatitis B, Polio, Measles, Herpes, Influenza, Rabies, Norwalk, Pox viruses, Enteroviruses, viral gastro-enteritis, etc, and over 400 strains of bacteria. Can be dissolved in water for soaking and wiping, or used in neat powder-form for decontaminating spills



Product No	Description
10285/01	Pack of 50 x 50g sachets
10285/02	Pack of 6 x 500g squeeze containers
10285/03	5kg drum

Azowipes

Product No	Description
10176/01	Non-woven wipers (220 x 210mm) impregnated with isopropyl alcohol for the safe disinfection of clinical and food areas. Pack of 200 wipers

Bench & Floor Protection

Kimwipes (7117)

3-ply steel blue wipers (508 x 457mm) provide clean absorbent covering for bench tops or use to mop up small spills. Dispensers follow



Product No	Qty
10160/01	Case of 12 rolls of 90 sheets

Dispenser (7049)

Durable silver-coloured steel roll holder for 10160/01 Kimwipe 7117 tissues

Product No	Description
10161/01	



Dispenser (7056)

White plastic-cased roll dispenser for 10160/01 Kimwipe 7117 tissues with tear edge for convenient one-handed dispensing

Product No	Description
10161/02	



Polytrap (295PE)

SCHLEICHER & SCHUELL highly absorbent paper sheets capable of soaking up 90ml/m² of water. Polyethylene-coated on one side to give optimum protection. Can be pre-soaked in disinfectant for use as an underlay for worktops where pathogens are present.

Product No	Description
12351/01	Pack of 50 sheets (480 x 600mm)
12351/02	50 metre roll, 400mm wide
12351/03	50 metre roll, 600mm wide
12351/04	50 metre roll, 1200mm wide

Polytrap (296PE)

SCHLEICHER & SCHUELL paper sheets similar to Polytrap 295PE, but with more than double the absorbency at 190ml/m² of water. Other sizes available on request.

Product No	Size
12351/10	Pack of 50 sheets (480 x 600mm)
12351/11	50 metre roll, 400mm wide
12351/12	50 metre roll, 600mm wide
12351/13	50 metre roll, 1200mm wide

Benchcote

WHATMAN absorbent paper sheets, laminated with polyethylene. Can be pre-soaked in disinfectant for use as an underlay.



Product No	Size
10074/01	Pack of 50 sheets (460 x 570mm)
10074/02	50 metre roll, 460mm wide
10074/03	50 metre roll, 920mm wide
10074/04	Pad of 50 sheets (460 x 570mm)

Benchcote Plus

WHATMAN absorbent paper sheets (500 x 600mm) similar to above, but capable of soaking up over 750ml/m²

Product No	Size
10075/01	Pack of 50 sheets
10075/02	50 metre roll, 600mm wide